

# Formstack Integration Data Guide



The Formstack Integration is supported for all Unlock, Medicom, and HealthAware HRAs. See additional tabs for any special variations that are HRA family specific (e.g. Unlock HRAs vs. Medicom Health HRAs, etc.). All data is written to the Lead object with a new Lead object created for each HRA completion.

## Lead Standard Fields

Lead Standard Fields are fields that are part of the standard Lead Object in Formstack. These fields will always populate if they are available and require no additional configuration. These fields are provided for all HRAs.

Formstack Lead Standard Field Name	Formstack Field Type	HRA Families	Unlock/HealthAware Field Name	Medicom Field Name	Expected Values	Description
Name	Name	All	first_name, last_name	first_name, last_name	1-100 character name	The name of the consumer. For HRAs that support caregivers answering on behalf of a person, it is NOT the name of the caregiver.
FirstName	Text	All	first_name	first_name	1-50 character first name	The first name of the consumer. For HRAs that support caregivers answering on behalf of a person, it is NOT the name of the caregiver.
LastName	Text	All	last_name	last_name	1-50 character last name	The last name of the consumer. For HRAs that support caregivers answering on behalf of a person, it is NOT the name of the caregiver.
Address	address	All	address_1 City Zip_code	address_2 State City State Zip_code	Full Address Name	The Full address of the consumer. Based on client customization, this field may be optional or hidden. If it is hidden, it will be empty.
street_address	Text	All	address_1	address_2 address_2	Address 1 & 2 Name	The street address of the consumer. Based on client customization, this field may be optional or hidden. If it is hidden, it will be empty.
City	Text	All	city	city	city name	The city of the consumer. Based on client customization, this field may be optional or hidden. If it is hidden, it will be empty.
state	Text	All	state	state	2 letter state initial in uppercase, ISO format ex: 'MN'	The 2 character state abbreviation of the US state for the consumer. Based on client customization, this field may be optional or hidden. The integration checks for StateCode first to determine if State and Country picklists have been enabled and if so, writes to the StateCode field, otherwise it adds the 2 letter state code into the state field directly.
zip_code	Text	All	zip_Code	zip_Code	5 digit string	The 5 digit zip code of the consumer. Based on client customization, this field may be optional or hidden. If it is hidden, it will be empty.
date_of_birth	datetime	All	date_of_birth	date_of_birth		The date of birth of the consumer. Format YYYYMMDD
phone	Phone	All	phone_1	phone_1	10 digit string	The 10 digit phone number of the consumer (US based number assumed).
utm_source	text	All	source	source	String from URL	This is the data that was passed to the HRA in the utm_source parameter. For Medicom HRAs, the s parameter is also able to be used.
Email	Email	All	email	email	email address	The email address of the consumer.
marketing_opt_in	checkbox	All	marketingOptIn	marketingOptIn	Unchecked = false Checked = true	Indicates if the consumer has opted-in to receive email marketing communications. If this option is not included in the HRA, the value is set based on parameters in the HRA Legal and Compliance section.
date_taken	Date	All	created_at	created_at	created	Date HRA was created. Format YYYYMMDD
Sex	Text	All	gender	gender	gender	The consumer's sex at birth. This should not be confused with gender. The Gender field names are still supported for continuity, but the data captured within the HRA is "Sex at Birth" and should not be used for any gender identity purposes.
client_data	Text	All	client_data	client_data	The "clientdata" parameter passed to the HRA as a URL parameter.	This is the "clientdata" parameter passed to the HRA as a URL parameter. Sometimes used as an external system identifier to link sessions to other systems.
ethnicity_ethnicity	Text	All	Ethnicity	Ethnicity	'WHITE' 'AFRICAN' 'HISPANIC' 'ASIAN_PACIFIC_ISLANDER' 'AMERICAN_INDIAN' 'OTHER' "	The ethnicity of the consumer from Medicom HRAs.
hra_type	Text	All	applicationTypeName	applicationTypeName	Acid Reflux Profiler Anxiety Profiler Back & Neck Pain Profiler Bladder Control Profiler Breast Cancer Risk Profiler	This is the type of HRA that was completed by the consumer.

- Breast Cancer Screening Assessment
- Colorectal Cancer Profiler
- Depression Profiler
- Diabetes Risk Profiler
- Healthy Weight-Range Profiler
- Heart Health Profiler
- Knee & Hip Profiler
- Lung Cancer Profiler
- PAD Symptoms Profiler
- Prostate Cancer Profiler
- Sleep Apnea Profiler
- Stroke Risk Profiler
- Weight-Loss Surgery Profiler
- Acid Reflux
- Anxiety
- Back & Neck
- Breast Cancer
- Bladder Control
- Colon Cancer
- COPD
- Depression
- Diabetes
- Heart
- Knee & Hip
- Lung Cancer
- Medical Weight Loss
- Peripheral Artery Disease
- Prostate Cancer
- Shoulder & Elbow
- Sleep
- Stroke
- Substance Use
- Weight Loss

bmi	Number	All	bmi	bmi	Decimal value representing BMI	This is the consumer's calculated BMI based on their answers.
height	Number	All	height	height	Integer value between 48 and 95 representing inches	This is the consumer's height in inches when they completed the HRA.
weight	Number	All	weight	weight	Integer value between 80 and 800 representing pounds	This is the consumer's weight in pounds when they completed the HRA.
has_pcp	Radio	Medicom	primaryCarePhysician_has	primaryCarePhysician_has	No Yes	The consumer's answer to "Do you have a (primary care name)?" where the label is customizable, but the default is primary care physician. Different from Medicom format because of a new "I don't know" option.
		Unlock HealthAware	primaryCarePhysician_has	primaryCarePhysician_has	Yes No I don't know	
pcp_last_visit	text	All	primaryCarePhysician_visited	primaryCarePhysician_visited		The consumer's answer to a question about the last time they visited their primary care physician for a checkup.
smoking_status	text	All	tobaccouse	tobaccouse	'NEVER' 'FORMER_MORE' 'FORMER_LESS' 'CURRENT' "	The smoking status of the consumer from Medicom HRAs.
utm_medium	text	All	medium	medium	String from URL	This is the data that was passed to the HRA in the utm_medium parameter.
utm_campaign	text	All	campaign	campaign	String from URL	This is the data that was passed to the HRA in the utm_campaign parameter.
utm_content	text	All	content	content	String from URL	This is the data that was passed to the HRA in the utm_content parameter.
utm_term	text	All	term	term	String from URL	This is the data that was passed to the HRA in the utm_term parameter.
phone_call_opt_in	Checkbox	All	phone_opt_in	phone_opt_in	Unchecked = false Checked = true	Indicates if the consumer has opted-in to receive a phone call about

						their results. If this option is not included in the HRA, the value is set based on parameters in the HRA Legal and Compliance section.
gender_pronouns	text	All	gender_pronouns	gender_pronouns		
hra_name	Text	All	app_type_name	app_type_name	Acid Reflux Profiler Anxiety Profiler Back & Neck Pain Profiler Bladder Control Profiler Breast Cancer Risk Profiler Breast Cancer Screening Assessment Colorectal Cancer Profiler Depression Profiler Diabetes Risk Profiler Healthy Weight-Range Profiler Heart Health Profiler Knee & Hip Profiler Lung Cancer Profiler PAD Symptoms Profiler Prostate Cancer Profiler Sleep Apnea Profiler Stroke Risk Profiler Weight-Loss Surgery Profiler Acid Reflux Anxiety Back & Neck Breast Cancer Bladder Control Colon Cancer COPD Depression Diabetes Heart Knee & Hip Lung Cancer Medical Weight Loss Peripheral Artery Disease Prostate Cancer Shoulder & Elbow Sleep Stroke Substance Use Weight Loss	This is the type of HRA that was completed by the consumer.
hra_location	Text	All	location	location	A string of the name of the HRA location label from the HRA Console.	This is the location name shown in the HRA Console.
Primary_Result	Text	All	risklevel	risklevel		Utilizes get_primary_result(@risk_level) to determine HRA active. Utilize Custom tabs to determine values for primary result by HRA Type
joint	text	All	joint	joint	'no_problems' 'mild_problems' 'moderate_problems' 'severe_problems'	
campaign_id	text	All	campaign_id	campaign_id		
custom_<custom question data name> (ex: custom_data_name)	HealthAware, Unlock	Unlock_Custom_Data_Name__c HealthAware_Custom_Data_Name__c	Depends on the custom question	Depends on the custom question	The custom question answer. You will need to add a field for each custom question you want to receive data for.	The custom question text as configured in the customization

# Formstack Integration Data Guide

The Formstack Integration is supported for all Unlock, Medcom, and HealthAware HRAs. See additional tabs for any special variations that are HRA family specific (e.g. Unlock HRAs vs. Medcom Health HRAs, etc.). All data is written to the Lead object with a new Lead object created for each HRA completion.



## Lead Custom Fields - Unlock HRA Specific

Lead Custom Fields are fields that can be added to the Lead Object in Formstack by an Administrator. These fields will be filled if they are available, but require additional set up by the client. Information is organized by HRA field name because we will look for multiple custom field names in order to store the data. Field names for parallel data types are aligned in order to easily migrate between HRA families.

HRA Data Name	Expected Data Values	Data Label	Description
primaryResult (from Acid Reflux HRA)	'SERIOUS_SYMPTOMS_REPORTED' 'DIAGNOSED_SYMPTOMATIC' 'UNDIAGNOSED_INCREASEDRISK' 'UNDIAGNOSED_SYMPTOMATIC' 'DIAGNOSED_ASYMPTOMATIC' 'UNDIAGNOSED_ASYMPTOMATIC'	Seek Medical Attention See a Doctor About Symptom Control See a Doctor for Evaluation Ask a Doctor About Your Symptoms Symptoms Well-Controlled GERD Not Likely	The primary result from the Acid Reflux HRA. It is backwards compatible with both Medcom and HealthAware HRAs.
primaryResult (from Anxiety HRA)	'RECOMMENDED' 'NOT_RECOMMENDED_SYMPTOMATIC' 'NOT_RECOMMENDED_ASYMPTOMATIC'	Anxiety assessment recommended Tell a doctor about your symptoms Keep a symptom journal	The primary result from the Anxiety HRA. It is backwards compatible with both Medcom and HealthAware HRAs.
backNeckPrimaryResult (from Back & Neck HRA)	'NO_IMPAIRMENT' 'MILD_IMPAIRMENT' 'MODERATE_IMPAIRMENT' 'SEVERE_IMPAIRMENT' 'OTHER_SYMPTOMS_REPORTED' 'INCONTINENCE_SYMPTOMS' 'MENINGITIS_OR_HEART_SYMPTOMS'	No Spine Impairment Mild Spine Impairment Moderate Spine Impairment Severe Spine Impairment See a Doctor Soon Seek Medical Attention Seek Immediate Medical Attention	The primary result from the Back & Neck HRA. It is backwards compatible with both Medcom and HealthAware HRAs.
primaryResult (from Bladder Control HRA)	'HAS_WORRISOME_URINARY_SYMPTOMS' 'GREATLY_BOTHERED_BY_STRESS_LEAKS' 'HAS_URGE_OR_STRESS_LEAKS_OR_FREQUENT_URINATION' 'HAS_LEAKS_NOT_RELATED_TO_URGE_OR_STRESS' 'NO_LEAKS_OR_FREQUENT_URINATION'	Prompt Medical Evaluation Recommended Treatable Control Problems Identified Treatable Bladder Problems Identified A Medical Evaluation Is Needed No Bladder Control Problems Reported	The primary result from the Bladder Control HRA. It is backwards compatible with both Medcom and HealthAware HRAs.
primaryResult (from Breast Cancer HRA)	'GENETIC_SCREENING_RECOMMENDED' 'GENETIC_RESCREENING_RECOMMENDED' 'SCREENING_RECOMMENDED' 'FOLLOW_UP_RECOMMENDED_SCREENING' 'DISCUSS_FUTURE_SCREENING' 'CONTINUE_ANNUAL_SCREENING' 'CONTINUE_RECOMMENDED_SCREENING'	Screening for Genetic Risk Factors Recommended Re-screening for Genetic Risk Factors Recommended Schedule Screening Now Follow Your Doctor's Recommended Screening Schedule Ask Your Doctor When You Should Begin Screening Continue Annual Screening Continue Recommended Screening Schedule	The primary result from the Breast Cancer HRA. It is backwards compatible with BreastCancerAware from HealthAware and the Breast Cancer Screening HRA from Medcom Health, but NOT the deprecated Breast Cancer Risk HRA from Medcom Health.
primaryResult (from Colon Cancer HRA)	'MONITORING_REC_EXISTING_CRC' 'DISCUSS_SYMPTOMS' 'DISCUSS_SCREENING_SCHEDULE' 'SCHEDULE_ROUTING_SCREENING_NOW' 'CONTINUE_ROUTINE_SCREENING' 'ABOVE_REC_SCREENING_AGE' 'BELOW_REC_SCREENING_AGE'	Personal Cancer History - Monitoring Recommended Discuss Symptoms with a Doctor Discuss Screening Schedule with a Doctor Schedule Routine Screening Now Continue Routine Screening Above Recommended Age for Screening Below Recommended Age for Screening	The primary result from the Colon Cancer HRA. It is backwards compatible with ColonCancerAware from HealthAware but NOT the Colon Cancer HRA from Medcom Health due to scientific improvements.
primaryResult (from COPD HRA)	'MANAGE_EXISTING_COPD' 'DISCUSS_EXISTING_COPD' 'SCREENING_RECOMMENDED' 'SCREENING_MAY_BE_RECOMMENDED' 'DISCUSS_SYMPTOMS'	Manage Your COPD With Your Doctor Talk to a Doctor About Your COPD Screening Recommended Screening May Be Recommended Tell a Doctor About Any Symptoms	The primary result from the COPD HRA. It is backwards compatible with the COPDAware HRA from HealthAware. Medcom Health does not have a compatible HRA.
primaryResult (from Depression HRA)	'RECOMMENDED' 'NOT_RECOMMENDED_SYMPTOMATIC' 'NOT_RECOMMENDED_ASYMPTOMATIC'	Depression Assessment Recommended Tell a Doctor About Your Symptoms Keep a Symptom Journal	The primary result from the Depression HRA. It is backwards compatible with both HealthAware and Medcom Health HRAs.
primaryResult (from Diabetes HRA)	'EXISTING_DIABETES' 'VERY_HIGH_RISK' 'HIGH_RISK' 'MODERATE_RISK' 'MORE_INFO_NEEDED' 'LOW_RISK'	Existing diabetes Very high risk High risk Moderate risk More info needed Low risk	The primary result from the Diabetes HRA. It is backwards compatible with DiabetesAware from HealthAware but NOT the Diabetes HRA from Medcom Health due to scientific improvements.
overallHeartRisk (from Heart HRA)	'EXISTING' 'HIGHNEARTERM' 'HIGHLONGTERM' 'MODERATE' 'INCREASED' 'LOW' 'UNKNOWN'	High Risk - Existing Condition High Near-term Risk High Long-term Risk Moderate Risk Increased Risk Low Risk Unknown Risk	The primary result from the Heart HRA. It is directly backwards compatible with HeartAware directly. A backwards-compatible option is also available for the Medcom Health Heart Health HRA data point for continuity. We strongly recommend upgrading to the new data type to benefit from additional stratification of the High Risk group.
medicomLegacyResult (from Heart HRA)	'UNKNOWN' 'LOW' 'MODERATE' 'HIGH' 'VERY_HIGH'	Unknown Risk Low Risk Moderate Risk High Risk High Risk	Converted data delivered for continuity purposes. Note that the data set includes 'VERY_HIGH' but is only displayed to the user as High Risk for user experience reasons. This is likely to be deprecated in the future, so clients should be migrating to the data type above.
jointAssessment (from Knee & Hip HRA)	'NO_PROBLEMS' 'MILD_PROBLEMS' 'MODERATE_PROBLEMS' 'SEVERE_PROBLEMS'	No Significant (joint) Problems Mild to Moderate (joint) Symptoms Moderate to Severe (joint) Symptoms Severe (joint) Symptoms	The primary result from the Knee & Hip HRA. It is backwards compatible with both HealthAware and Medcom Health HRAs.
primaryResult (from Lung Cancer HRA)	'RECOMMENDED' 'NOT_RECOMMENDED_SMOKER' 'DISCUSS' 'NOT_RECOMMENDED_QUIT' 'NOT_RECOMMENDED_YOUNG' 'NOT_RECOMMENDED_OLD' 'NOT_RECOMMENDED_NEVER'	Screening May Be Recommended Screening Not Currently Recommended Discuss Screening with a Doctor Former Smoker -- Screening Not Recommended Discuss Screening at Age 50 Screening Not Recommended Due to Age Never Smoked -- Screening Not Recommended	The primary result from the Lung Cancer HRA. It is backwards compatible with LungCancerAware from HealthAware but NOT the Lung Cancer HRA from Medcom Health due to scientific improvements.
primaryResult (from Medical Weight Loss HRA)	'LIKELY_ELIGIBLE' 'POSSIBLY_ELIGIBLE' 'NOT_ELIGIBLE'	Likely Eligible Possibly Eligible Not Eligible	The primary result from the Medical Weight Loss HRA.
screeningRecommendation (from Peripheral Artery Disease HRA)	'RECOMMENDED_SYMPTOMATIC' 'RECOMMENDED' 'INFO_NEEDED' 'NOT_RECOMMENDED'	Seek Medical Attention Soon Screening Recommended More Information Needed No Screening Recommended	The primary result from the Peripheral Artery Disease HRA. It is backwards compatible with both HealthAware and Medcom Health HRAs.
primaryResult (from Prostate Cancer HRA)	'EXISTING_PCA' 'HAS_EXISTING_SCREENING_PLAN' 'DISCUSS_SCREENING' 'DISCUSS_SCREENING_AT_AGE_FORTY' 'DISCUSS_SCREENING_AT_AGE_FORTYFIVE'	Follow Your Doctor's Recommendation Follow Your Existing Plan Discuss Screening Options Discuss Screening Options at Age 40 Discuss Screening Options at Age 45	The primary result from the Prostate Cancer HRA. It is backwards compatible with both HealthAware and Medcom Health HRAs.

primaryResult (from Shoulder & Elbow HRA)	'DISCUSS_SCREENING_AT_AGE_FIFTY' 'NO_SYMPTOMS' 'MILD_SYMPTOMS' 'MODERATE_SYMPTOMS' 'SEVERE_SYMPTOMS'	Discuss Screening Options at Age 50 No Significant (joint) Symptoms Mild to Moderate (joint) Symptoms Moderate to Severe (joint) Symptoms Severe (joint) Symptoms	The primary result from the Unlock Shoulder & Elbow HRA.
primaryResult (from Sleep HRA)	'HIGH_RISK' 'LOW_RISK' 'UNKNOWN_RISK'	High Sleep Apnea Risk Low Sleep Apnea Risk Unknown Risk for Sleep Apnea	The primary result from the Sleep HRA. It is backwards compatible with both HealthAware and Medicom Health HRAs.
primaryResult (from Stroke HRA)	'HIGH2' 'HIGH1' 'MODERATE' 'LOW'	High Risk - Existing Condition High Risk Moderate Risk Low Risk	The primary result from the Stroke HRA. It is backwards compatible with HealthAware. For Medicom Health HRA, the data types are also the same, but the label used within the HRA for 'HIGH2' was previously Very High Risk. This was changed HealthAware and Unlock HRAs to clarify that it is tied to an existing condition.
primaryResult (from Substance Use HRA)	'SEEK_IMMEDIATE_HELP' 'EVALUATION_RECOMMENDED' 'LEARN_ABOUT_TREATMENT' 'NOTIFY_PCP_OF_FUTURE_CHANGES'	Get Confidential Help Today See a Health Care Professional Learn More About Treatment Options Tell a Health Care Provider If Things Change	The primary result from the Substance Use HRA.
primaryResult (from Weight HRA)	'POSSIBLE_SURGERY_CANDIDATE' 'SCREENING_NEEDED' 'WEIGHT_LOSS_RECOMMENDED' 'RECOMMEND_QOL_FOLLOWUP' 'WEIGHT_LOSS_MAY_BE_RECOMMENDED' 'WEIGHT_LOSS_NOT_RECOMMENDED'	Possible Surgery Candidate Screening Needed to Determine Eligibility Weight Loss Recommended Recommend "Quality of Life" Follow-up Weight Loss May Be Recommended Weight Loss Not Recommended	The primary result from the Weight HRA. It is directly backwards compatible with HeartAware directly. A backwards-compatible option is also available for the Medicom Health Weight-Loss Surgery HRA data point for continuity. We strongly recommend upgrading to the new data type to benefit from additional stratification of the data.
medicomLegacyResult (from Weight HRA)	'MORE_INFO' 'NO' 'YES'	More information needed to determine a recommendation You don't appear to be a candidate for surgery. You may be a candidate for weight-loss surgery.	Converted data delivered for continuity purposes. This is likely to be deprecated in the future, so clients should be migrating to the data type above.

# Formstack Integration Data Guide



The Formstack Integration is supported for all Unlock, Medicom, and HealthAware HRAs. See additional tabs for any special variations that are HRA family specific (e.g. Unlock HRAs vs. Medicom Health HRAs, etc.). All data is written to the Lead object with a new Lead object created for each HRA completion.

## Lead Custom Fields - Medicom Health Specific

Lead Custom Fields are fields that can be added to the Lead Object in Formstack by an Administrator. These fields will be filled if they are available, but require additional set up by the client. Information is organized by HRA field name because we will look for multiple custom field names in order to store the data. Field names for parallel data types are aligned in order to easily migrate between HRA families.

HRA Data Name	Expected Data Values	Data Label	Description/HRA Type
symptomAssessment	'DIAGNOSED_SYMPTOMATIC'	See a Doctor About Symptom Control	Acid Reflux
	'DIAGNOSED_ASYMPTOMATIC'	Symptoms Well-Controlled	
	'UNDIAGNOSED_INCREASEDRISK'	See a Doctor for Evaluation	
	'UNDIAGNOSED_SYMPTOMATIC'	Ask a Doctor About Your Symptoms	
	'UNDIAGNOSED_ASYMPTOMATIC'	Acid Reflux Disease Not Likely	
primaryResult	'RECOMMENDED'	Anxiety assessment recommended	Anxiety
	'NOT_RECOMMENDED_SYMPTOMATIC'	Tell a doctor about your symptoms	
	'NOT_RECOMMENDED_ASYMPTOMATIC'	Keep a symptom journal	
backNeckPrimaryResult	'NO_IMPAIRMENT'	No Spine Impairment	Back & Neck
	'MILD_IMPAIRMENT'	Mild Spine Impairment	
	'MODERATE_IMPAIRMENT'	Moderate Spine Impairment	
	'SEVERE_IMPAIRMENT'	Severe Spine Impairment	
	'OTHER_SYMPTOMS_REPORTED'	See a Doctor Soon	
	'INCONTINENCE_SYMPTOMS'	Seek Medical Attention	
	'MENINGITIS_OR_HEART_SYMPTOMS'	Seek Immediate Medical Attention	
bladderControlAssessment	'HAS_WORRISOME_URINARY_SYMPTOMS'	Prompt Medical Evaluation Recommended	Bladder Control
	'GREATLY_BOTHERED_BY_STRESS_LEAKS'	Treatable Control Problems Identified	
	'HAS_URGE_OR_STRESS_LEAKS_OR_FREQUENT_URINATION'	Treatable Bladder Problems Identified	
	'HAS_LEAKS_NOT_RELATED_TO_URGE_OR_STRESS'	A Medical Evaluation Is Needed	
	'NO_LEAKS_OR_FREQUENT_URINATION'	No Bladder Control Problems Reported	
breastCancer5YearRiskCategory	'HIGH'	High 5-Year Risk = xx.xx% (>=3%)	Breast Cancer Risk
	INCREASED'	Increased 5-Year Risk = x.xx% (>=1.66% AND < 3%)	
	'LOW'	Low 5-Year Risk = x.xx% (<1.66%)	
	'UNKNOWN_AGE'	Unknown 5-Year Risk - Cannot Be Calculated (due to age)	
	'UNKNOWN_HISTORY'	Unknown 5-Year Risk - Incomplete Family History	
primaryResult	'GENETIC_SCREENING_RECOMMENDED'	Screening for Genetic Risk Factors Recommended	Breast Cancer Screening
	'GENETIC_RESCREENING_RECOMMENDED'	Re-screening for Genetic Risk Factors Recommended	
	'SCREENING_RECOMMENDED'	Schedule Screening Now	
	'FOLLOW_RECOMMENDED_SCREENING'	Follow Your Doctor's Recommended Screening Schedule	
	'DISCUSS_FUTURE_SCREENING'	Ask Your Doctor When You Should Begin Screening	
	'CONTINUE_ANNUAL_SCREENING'	Continue Annual Screening	
	'CONTINUE_RECOMMENDED_SCREENING'	Continue Recommended Screening Schedule	
screening	'FOLLOW_DOC_REC'	Follow Screening Schedule Set by Doctor	Colorectal Cancer
	'REVIEW_W_DOC'	Review Screening Schedule with a Doctor	
	'RECOMMENDED'	Screening Recommended	
	'START_AT_40'	Start Screening at Age 40	
	'START_AT_45'	Start Screening at Age 45	
	'DISCUSS_SCHEDULE'	Discuss Screening Schedule with a Doctor	
	'RECOMMENDED'	Depression Assessment Recommended	
'NOT_RECOMMENDED_SYMPTOMATIC'	Tell a Doctor About Your Symptoms		
'NOT_RECOMMENDED_ASYMPTOMATIC'	Keep a Symptom Journal		
diabetesRisk	'NA'	Not Applicable	Diabetes
	'LOW'	Low Risk	
	'MODERATE'	Moderate Risk	
	'HIGH'	High Risk	
	'VERY_HIGH'	Very High Risk	
weightProfile	'UNEXPLAINED_WEIGHT_LOSS'	See a Doctor About Your Unexplained Weight Loss	Healthy Weight Range
	'OBESE_LOSE_WEIGHT'	Current: Obesity; Healthy: Below ### lbs.	
	'OVERWEIGHT_LOSE_WEIGHT'	Current: Overweight; Healthy: Below ### lbs.	

	'OVERWEIGHT_MAINTAIN_WEIGHT'	Currently in the Overweight Range	
	'NORMAL_MAINTAIN_WEIGHT'	Currently in the Healthy Weight Range	
	'UNDERWEIGHT_GAIN_WEIGHT'	Currently in the Underweight Range	
riskCategory	'UNKNOWN'	Unknown Risk	Heart Health
	'LOW'	Low Risk	
	'MODERATE'	Moderate Risk	
	'HIGH'	High Risk	
	'VERY_HIGH'	High Risk (<- This is correct)	
jointAssessment	'NO_PROBLEMS'	No Significant {joint} Problems	Knee & Hip
	'MILD_PROBLEMS'	Mild to Moderate {joint} Symptoms	
	'MODERATE_PROBLEMS'	Moderate to Severe {joint} Symptoms	
	'SEVERE_PROBLEMS'	Severe {joint} Symptoms	
screening	'RECOMMENDED'	Screening May Be Recommended	Lung Cancer
	'NOT_RECOMMENDED_SMOKER'	Not Currently Recommended	
	'DISCUSS'	Discuss Screening Options with a Doctor	
	'NOT_RECOMMENDED_QUIT'	Former Smoker -- Not Recommended	
	'NOT_RECOMMENDED_YOUNG'	Discuss Screening Options at Age 55	
	'NOT_RECOMMENDED_OLD'	Not Recommended Due to Age	
	'NOT_RECOMMENDED_NEVER'	Never Smoked -- Not Recommended	
screeningRecommendation	'RECOMMENDED_SYMPTOMATIC'	Seek Medical Attention Soon	PAD
	'RECOMMENDED'	Screening Recommended	
	'INFO_NEEDED'	More Information Needed	
	'NOT_RECOMMENDED'	No Screening Recommended	
primaryResult	'EXISTING_PCA'	Follow Your Doctor's Recommendation	Prostate Cancer
	'HAS_EXISTING_SCREENING_PLAN'	Follow Your Existing Plan	
	'DISCUSS_SCREENING'	Discuss Screening Options	
	'DISCUSS_SCREENING_AT_AGE_FORTY'	Discuss Screening Options at Age 40	
	'DISCUSS_SCREENING_AT_AGE_FORTYFIVE'	Discuss Screening Options at Age 45	
	'DISCUSS_SCREENING_AT_AGE_FIFTY'	Discuss Screening Options at Age 50	
sleepApneaRiskCategory	'HIGH_RISK'	High Sleep Apnea Risk	Sleep Apnea
	'LOW_RISK'	Low Sleep Apnea Risk	
	'UNKNOWN_RISK'	Unknown Risk for Sleep Apnea	
strokeRiskCategory	'HIGH2'	Very High Risk	Stroke
	'HIGH1'	High Risk	
	'MODERATE'	Moderate Risk	
	'LOW'	Low Risk	
weightLossSurgeryRecommendation	'MORE_INFO'	More information needed to determine a recommendation.	Weight Loss Surgery
	'NO'	You don't appear to be a candidate for surgery.	
	'YES'	You may be a candidate for weight-loss surgery.	

# Formstack Integration Data Guide

The Formstack Integration is supported for all Unlock, Medicom, and HealthAware HRAs. See additional tabs for any special variations that are HRA family specific (e.g. Unlock HRAs vs. Medicom Health HRAs, etc.). All data is written to the Lead object with a new Lead object created for each HRA completion. For more information about this integration, please refer to the "SF Health Cloud - Unlock HRAs Integration" document.



## Lead Custom Fields - HealthAware Specific

Lead Custom Fields are fields that can be added to the Lead Object in Formstack by an Administrator. These fields will be filled if they are available, but require additional set up by the client. Information is organized by HRA field name because we will look for multiple custom field names in order to store the data. Field names for parallel data types are aligned in order to easily migrate between HRA families.

HRA Data Name	Expected Data Values	Data Label	Description/HRA Type
primaryResult	'SERIOUS_SYMPTOMS_REPORTED'	Seek Medical Attention	GERDAware
	'DIAGNOSED_SYMPTOMATIC'	See a Doctor About Symptom Control	
	'UNDIAGNOSED_INCREASEDRISK'	See a Doctor for Evaluation	
	'UNDIAGNOSED_SYMPTOMATIC'	Ask a Doctor About Your Symptoms	
	'DIAGNOSED_ASYMPTOMATIC'	Symptoms Well-Controlled	
primaryResult	'UNDIAGNOSED_ASYMPTOMATIC'	GERD Not Likely	AnxietyAware
	'RECOMMENDED'	Anxiety assessment recommended	
PrimaryResult	'NOT_RECOMMENDED_SYMPTOMATIC'	Tell a doctor about your symptoms	SpineAware
	'NOT_RECOMMENDED_ASYMPTOMATIC'	Keep a symptom journal	
PrimaryResult	'NO_IMPAIRMENT'	No Spine Impairment	SpineAware
	'MILD_IMPAIRMENT'	Mild Spine Impairment	
	'MODERATE_IMPAIRMENT'	Moderate Spine Impairment	
	'SEVERE_IMPAIRMENT'	Severe Spine Impairment	
	'OTHER_SYMPTOMS_REPORTED'	See a Doctor Soon	
	'INCONTINENCE_SYMPTOMS'	Seek Medical Attention	
primaryResult	'MENINGITIS_OR_HEART_SYMPTOMS'	Seek Immediate Medical Attention	BladderAware
	'HAS_WORRISOME_URINARY_SYMPTOMS'	Prompt Medical Evaluation Recommended	
	'GREATLY_BOTHERED_BY_STRESS_LEAKS'	Treatable Control Problems Identified	
	'HAS_URGE_OR_STRESS_LEAKS_OR_FREQUENT_URINATION'	Treatable Bladder Problems Identified	
	'HAS_LEAKS_NOT_RELATED_TO_URGE_OR_STRESS'	A Medical Evaluation Is Needed	
primaryResult	'NO_LEAKS_OR_FREQUENT_URINATION'	No Bladder Control Problems Reported	BreastCancerAware
	'GENETIC_SCREENING_RECOMMENDED'	Screening for Genetic Risk Factors Recommended	
	'GENETIC_RESCREENING_RECOMMENDED'	Re-screening for Genetic Risk Factors Recommended	
	'SCREENING_RECOMMENDED'	Schedule Screening Now	
	'FOLLOW_RECOMMENDED_SCREENING'	Follow Your Doctor's Recommended Screening Schedule	
	'DISCUSS_FUTURE_SCREENING'	Ask Your Doctor When You Should Begin Screening	
	'CONTINUE_ANNUAL_SCREENING'	Continue Annual Screening	
'CONTINUE_RECOMMENDED_SCREENING'	Continue Recommended Screening Schedule		
primaryResult	'MONITORING_REC_EXISTING_CRC'	Personal Cancer History - Monitoring Recommended	ColonCancerAware
	'DISCUSS_SYMPTOMS'	Discuss Symptoms with a Doctor	
	'DISCUSS_SCREENING_SCHEDULE'	Discuss Screening Schedule with a Doctor	
	'SCHEDULE_ROUTING_SCREENING_NOW'	Schedule Routine Screening Now	
	'CONTINUE_ROUTINE_SCREENING'	Continue Routine Screening	
	'ABOVE_REC_SCREENING_AGE'	Above Recommended Age for Screening	
primaryResult	'BELOW_REC_SCREENING_AGE'	Below Recommended Age for Screening	COPDAware
	'MANAGE_EXISTING_COPD'	Manage Your COPD With Your Doctor	
	'DISCUSS_EXISTING_COPD'	Talk to a Doctor About Your COPD	
	'SCREENING_RECOMMENDED'	Screening Recommended	

primaryResult	'SCREENING_MAY_BE_RECOMMENDED'	Screening May Be Recommended	
	'DISCUSS_SYMPTOMS'	Tell a Doctor About Any Symptoms	
	'RECOMMENDED'	Depression Assessment Recommended	DepressionAware
	'NOT_RECOMMENDED_SYMPTOMATIC'	Tell a Doctor About Your Symptoms	
primaryResult	'NOT_RECOMMENDED_ASYMPTOMATIC'	Keep a Symptom Journal	
	'EXISTING_DIABETES'	Existing diabetes	DiabetesAware
	'VERY_HIGH_RISK'	Very high risk	
	'HIGH_RISK'	High risk	
	'MODERATE_RISK'	Moderate risk	
	'MORE_INFO_NEEDED'	More info needed	
overallHeartRisk	'LOW_RISK'	Low risk	
	'EXISTING'	High Risk - Existing Condition	HeartAware
	'HIGHNEARTERM'	High Near-term Risk	
	'HIGHLONGTERM'	High Long-term Risk	
	'MODERATE'	Moderate Risk	
	'INCREASED'	Increased Risk	
jointAssessment	'LOW'	Low Risk	
	'UNKNOWN'	Unknown Risk	
	'NO_PROBLEMS'	No Significant {joint} Problems	JointAware - Knee & Hip
	'MILD_PROBLEMS'	Mild to Moderate {joint} Symptoms	
	'MODERATE_PROBLEMS'	Moderate to Severe {joint} Symptoms	
primaryResult	'SEVERE_PROBLEMS'	Severe {joint} Symptoms	
	'RECOMMENDED'	Screening May Be Recommended	LungCancerAware
	'NOT_RECOMMENDED_SMOKER'	Screening Not Currently Recommended	
	'DISCUSS'	Discuss Screening with a Doctor	
	'NOT_RECOMMENDED_QUIT'	Former Smoker -- Screening Not Recommended	
	'NOT_RECOMMENDED_YOUNG'	Discuss Screening at Age 50	
screeningRecommendation	'NOT_RECOMMENDED_OLD'	Screening Not Recommended Due to Age	
	'NOT_RECOMMENDED_NEVER'	Never Smoked -- Screening Not Recommended	
	'RECOMMENDED_SYMPTOMATIC'	Seek Medical Attention Soon	VascularAware
	'RECOMMENDED'	Screening Recommended	
	'INFO_NEEDED'	More Information Needed	
	'NOT_RECOMMENDED'	No Screening Recommended	
primaryResult	'EXISTING_PCA'	Follow Your Doctor's Recommendation	ProstateCancerAware
	'HAS_EXISTING_SCREENING_PLAN'	Follow Your Existing Plan	
	'DISCUSS_SCREENING'	Discuss Screening Options	
	'DISCUSS_SCREENING_AT_AGE_FORTY'	Discuss Screening Options at Age 40	
	'DISCUSS_SCREENING_AT_AGE_FORTYFIVE'	Discuss Screening Options at Age 45	
primaryResult	'DISCUSS_SCREENING_AT_AGE_FIFTY'	Discuss Screening Options at Age 50	
	'HIGH_RISK'	High Sleep Apnea Risk	SleepAware
	'LOW_RISK'	Low Sleep Apnea Risk	
primaryResult	'UNKNOWN_RISK'	Unknown Risk for Sleep Apnea	
	'HIGH2'	High Risk - Existing Condition	StrokeAware
	'HIGH1'	High Risk	
	'MODERATE'	Moderate Risk	
primaryResult	'LOW'	Low Risk	
	'POSSIBLE_SURGERY_CANDIDATE'	Possible Surgery Candidate	WeightAware
	'SCREENING_NEEDED'	Screening Needed to Determine Eligibility	
	'WEIGHT_LOSS_RECOMMENDED'	Weight Loss Recommended	
	'RECOMMEND_QOL_FOLLOWUP'	Recommend "Quality of Life" Follow-up	
	'WEIGHT_LOSS_MAY_BE_RECOMMENDED'	Weight Loss May Be Recommended	

'WEIGHT\_LOSS\_NOT\_RECOMMENDED'

Weight Loss Not Recommended